



Product Specification Sheet

EpiSep Urine Conditioner®

FOR RESEARCH USE ONLY: This product is not for diagnostic or therapeutic use.

CATALOG #: R2125-1

DESCRIPTION: Buffered pH 6.5, urine pre-treatment that reduces crystal formation and breaks up chylomicrons that are commonly found in urine.

SUPPLIED AS: 100mL, ready to use liquid.

STABILITY: Stable when stored at 18-24°C. Recommended storage upon receipt is 18-24°C. **KEEP CONTAINER CLOSED. KEEP AWAY FROM HEAT, SPARKS, AND OPEN FLAME.**

APPLICATION SHEETS:

EpiSep Preparation and Isolation of Urothelial Cells from Preserved Urine for use with FISH. (Ref. Doc#ML080612-1)

MATERIAL SAFETY DATA:

When handling this material Standard Laboratory Practices or Proper Protective Equipment should be followed. Avoid breathing vapor. Avoid contact with eyes, skin and clothing. Contains Sodium Azide as a preservative. Although, the quantity of sodium azide (0.09%) is very small, measures should be taken to avoid skin and eye contact, inhalation and ingestion. Sodium Azide (NaN₃) may react with lead and copper plumbing to form potentially explosive metal oxides. Upon disposal, flush with a large volume of water to prevent Azide build-up.

